

INSTRUCTIONS FOR USE

en

I. COMPOSITION

Cross-linked hyaluronic acid	20 mg/mL
Lidocaine hydrochloride	3 mg/mL
Phosphate buffered saline pH 7	qs ad 1 mL

II. DESCRIPTION

Restylane® Kysse™ is a sterile, biodegradable, transparent gel of non-animal cross-linked hyaluronic acid with the addition of lidocaine hydrochloride 3 mg/mL. The gel is supplied in a prefilled plastic syringe. The contents of the syringe are sterilized using moist heat. The syringe is packaged individually in a blister, with two 30G x 1/2" Ultra thin wall needles. The needles have been sterilized using irradiation. The product is for single use only. To ensure traceability the package includes patient record labels that should be attached to patient records.

III. INTENDED USE

The product is intended to restore or augment the volume of the lips. It should be injected into the submucosal layer of the lip. Lidocaine is added to the formulation to diminish the pain resulting from the injection during the treatment.

The product is only intended to be used by authorized personnel in accordance with local legislation, trained in the appropriate injection techniques. Before the first treatment session, it is recommended to contact your local Galderma representative or Restylane distributor for information about training opportunities.

IV. MODE OF ACTION

The product adds volume to the tissue thereby enhancing the lips to the desired level. The volume and lifting capacity originate from the ability of cross-linked hyaluronic acid to bind water.

V. CONTRAINDICATIONS

- Patients presenting with known allergy to hyaluronic acid filler or amide local anaesthetics.
- Patients presenting with porphyria.

VI. WARNINGS

- Do not use where there is active disease, such as inflammation, infection or tumours, in or near the intended treatment site.
- Do not inject intravascularly. As for other injectable medical devices, inadvertent injection into or next to blood vessels could potentially lead to vascular occlusion or compression, ischemia and necrosis.
- Do not use in patients with bleeding disorders or in patients who are taking thrombolytics or anticoagulants.
- Do not inject this product into an area where an implant other than hyaluronic acid has been placed.
- Do not resterilize.
- Do not mix with other products.

VII. PRECAUTIONS FOR USE

General considerations relevant to injectable medical devices

- Knowledge of the anatomy of treatment site and special caution are required in order to avoid perforation or compression of vessels, nerves and other vulnerable structures.
- Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be observed.
- Patients using immunosuppressants are not suitable candidates for treatment.
- Special caution should be exercised in treating patients with a tendency to form hypertrophic scars or any other healing disorders.
- Injection procedures can lead to reactivation of latent or subclinical herpes viral infections.
- Patients who are using substances that affect platelet function, such as aspirin and nonsteroidal anti-inflammatory drugs may, as with any injection, experience increased bruising or bleeding at injection sites.
- Patients with unattainable expectations are not suitable candidates for treatment.

- Do not use the product if the package is damaged.
- Do not use the product if the contents of the syringe are cloudy.

Specific considerations relevant to the use of this product

- Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High doses of lidocaine (more than 400 mg) can cause acute toxic reactions manifesting as symptoms affecting the central nervous system and cardiac conduction.
- Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type local anaesthetics e.g., certain anti-arrhythmics, since the systemic toxic effects can be additive.
- Lidocaine should be used cautiously in patients with epilepsy, impaired cardiac conduction, severely impaired hepatic function or severe renal dysfunction.
- If the product is injected too superficially this may result in visible lumps and/or bluish discoloration.
- The patient must avoid exposing the treated area to heat (sun bathing, sauna, steam baths, etc.) or extreme cold until any signs of local inflammation have disappeared.
- If laser treatment, chemical peeling or any other procedure based on active dermal response is performed after treatment with this product there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if the product is administered before the skin has healed completely after such a procedure.
- At each treatment session a maximum dosage of 1.5 mL per treatment site is recommended.
- This product has not been tested in pregnant or breastfeeding women.
- Do not use in children.

VIII. ADVERSE EVENTS

Patients must be informed of the potential risks and adverse events related to the injection procedure and to the use of this product.

The following post marketing adverse events have been reported (non-exhaustive list): angioedema, atrophy/scarring, blisters, bruising, erythema, hypersensitivity, infection, mass, neurological symptoms such as paraesthesia, pain/tenderness, papules/nodules, pruritus, reactivation of herpes infection, short duration of effect, and swelling.

Other potential adverse events, that have been reported following treatment of lips with hyaluronic acid gels in general and may occur when using the product include the following: abscess, acne, capillary disorders such as telangiectasia, dermatitis, device dislocation, discoloration, granuloma, induration, inflammation, and ischemia/necrosis.

Injection related adverse events such as bruising, erythema, itching, swelling, pain and tenderness generally resolve spontaneously within one week after injection.

Vascular compromise may occur due to an inadvertent intravascular injection or as a result of vascular compression associated with implantation of any injectable product. This may manifest as ischemia or necrosis at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolisation.

Symptoms of inflammation at the implant site commencing either shortly after injection or after a delay of up to several weeks have been reported. In case of unexplained inflammatory reactions infections should be excluded and treated if necessary since inadequately treated infections may progress into complications such as abscess formation. Treatment using only oral corticosteroids without concurrent antibiotic treatment is not recommended. In case of persistent or recurrent inflammatory symptoms, consider removal of the product by aspiration/drainage, extrusion or enzymatic degradation (use of hyaluronidase has been described in scientific publications). Before any removal procedure is performed, the swelling may be reduced by using e.g. NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product.

For patients who have experienced medically important adverse events, a decision for retreatment should take into consideration the cause and severity of previous reactions.

For reporting adverse events contact your local Galderma representative or distributor for this product.

IX. TREATMENT PROCEDURE

Inform the patient about the precautions to be taken, the expected result and the possible adverse events.

It is important to use a sterile, appropriate needle or blunt cannula. Suitable needles (30G x 1/2" Ultra thin wall) are supplied with the syringe in the blister pack. As an alternative, a blunt thin walled cannula with a recommended size of 25-27G can be used.

Assembling the needle/cannula and syringe:

- use surgical gloves
- unscrew the protective cap from the stopper
- carefully remove the stopper from the syringe
- firmly screw the needle/cannula with its shield onto the tip of the syringe
- remove the shield just before injection

Cleanse the area to be treated with an antiseptic and allow it to dry before injection.

To avoid breakage of the needle/cannula, do not attempt to bend or otherwise manipulate it before or during treatment.

Before injecting the product, depress the plunger rod carefully until a small droplet is visible at the tip of the needle/cannula.

Align the bevel of the needle by turning the syringe on its axis. If a blunt cannula is used, an entry point is made in the skin, for example with a sharp needle of appropriate size.

Aspiration is recommended prior to injection in order to reduce the risk of inadvertent injection into a blood vessel. Inject slowly. During injection, keep the side hole of the cannula facing downwards, away from the skin surface. Inject the gel by gently pressing down on the plunger rod with the thumb or palm of the hand.

Choose from a variety of injection techniques, i.e. serial puncture, linear threading or cross-hatching.

It is recommended to change needle/cannula for each new treatment site. Defects should be fully corrected, but not overcorrected, at each treatment session.

If "blanching" of the skin is observed as a result of over-superficial injection, the whitened area should be massaged gently until it returns to a normal colour.

Gently massage the treated area after injection.

The syringes and needles/cannulas must be discarded immediately after use in accordance with accepted medical practice and applicable national, local or institutional guidelines. The product shall not be reused due to risk for contamination of unused material and the associated risk including infection. Standard precautions apply when handling the needles. Needles should be disposed in a container dedicated for sharp devices.

X. SHELF LIFE AND STORAGE



Do not use after the expiry date indicated on package. Store up to 25°C. Protect from freezing and sunlight.

XI. MANUFACTURER

Q-Med AB, Seminariegatan 21, SE-75228 Uppsala, Sweden
Phone +46(0)18 474 90 00, Fax +46(0)18 474 90 01
www.q-med.com, e-mail: info.q-med@galderma.com

Restylane, Emervel and Galderma are trademarks of Nestlé Skin Health S.A.

Symbols on the packaging

		
en	CE marking in compliance with directive 93/42/EEC concerning medical devices. The accompanied identification number is the number of the Notified Body.	Do not use if package is damaged.
hr	CE-oznaka uskladenost s odredbama Smjernice 93/42/EEC koja se odnosi na medicinske uređaje. Popratni identifikacijski broj je broj prijavljenog tijela.	Ne upotrebljavajte ako je pakiranje oštećeno.
cs	Označení CE dle směrnice 93/42/EHS o zdravotnických prostředcích. Identifikační číslo vedle značky je registračním číslem certifikující notifikované osoby.	Nepoužívejte, je-li obal poškozen.
no da	CE-merket i henhold til EU-direktiv 93/42/EØF om medisinsk utstyr. Tilhørende identifikasjonsnummer er nummer på meldt/bemyndiget organ.	Må ikke benyttes dersom emballasjen er skadet.
nl	CE-markering in overeenstemming met Richtlijn 93/42/EEC met betrekking tot medische hulpmiddelen. Het bijbehorende identificatienummer is het nummer van de aangewezen instantie.	Niet gebruiken indien de verpakking beschadigd is.
fi	Lääkintälaite-direktiivin 93/42/ETY mukainen CE-merkintä. Merkinän numero on ilmoitetun laitoksen numero.	Älä käytä, jos pakkaus on vahingoittunut.
fr	Marquage CE conforme à la directive 93/42/CEE sur les dispositifs médicaux. Le numéro d'identification est le numéro de l'organisme notifié.	Ne pas utiliser si la boîte est endommagée.
de	CE-Kennzeichnung gemäß Richtlinie 93/42/EWG über Medizinprodukte. Die angegebene Identifikationsnummer ist die Nummer der benannten Stelle.	Nicht verwenden, wenn die Verpackung Schäden aufweist.
sr	CE oznaka u skladu sa Direktivom 93/42/EEC o medicinskim uređajima. Prateći identifikacioni broj je broj Ovlašćenog tela.	Ne koristiti ako je pakovanje oštećeno.
sk	CE označujúca zhodu so Smernicou 93/42/EEC o zdravotníckych pomôckach. Správádzajúce identifikačné číslo je číslo autorizovaného orgánu.	Nepoužívajte, ak je balenie poškodené.
sv	CE märkt enligt direktivet 93/42/EEG om medicintekniska produkter. Identifikationsnumret är numret på anmält organ.	Får inte användas om förpackningen är skadad.